


How to use long reads for de-novo assembly


Wolfram Höps
Radboudumc Nijmegen

Complex SVs can be hard to interpret from aligned reads




Several Structural Variants are visible, but which genes are affected?

Complex loci can be hard to interpret from aligned reads



A multi-Mbp-sized deletion is visible, but where are the breakpoints?

Alignment




Read alignment != de-novo assembly

Alignment


≠

De-novo assembly




Read alignment != de-novo assembly

Alignment




≠

De-novo assembly




Read alignment != de-novo assembly



Requires Reference Genome
Short variant calling in non-complex regions
Computationally light (hours)
Automated, easy to use

No Reference required (though sometimes useful)
Long, complex variants in hard-to-map regions
Computationally heavy (hours - days)
Pitfalls, hard to use

The telomere-to-telomere reference **assembly**



red:
regions
previously
(hg38)
unresolved

Sequencing all 70,000
eukaryotic species of
Britain and Ireland




[Portal.darwintreeoflife.org](https://portal.darwintreeoflife.org)

Sequencing all 70,000 eukaryotic species of Britain and Ireland



[Portal.darwintreeoflife.org](https://portal.darwintreeoflife.org)



Human Pangenome Reference Consortium:
232 individuals fully assembled (-> 464 'genomes')

The assembly 'hype' is technology-driven




The assembly 'hype' is technology-driven

Short reads

DNA molecule



Short Reads




Puzzling is very hard - and complex or repetitive regions can barely be resolved

Long reads

DNA molecule




Long Reads



Puzzling becomes much easier

De-novo assembly of whole human genomes or smaller loci is now feasible for well-equipped labs



Typical applications in human genetics:

- resolve structural variations (SVs) longer than the read length
- resolve variants in or near unmappable or highly diverse genomic regions
- resolve complex genomic regions

Required

>15-fold accurate long reads per haplotype is needed.

- PacBio HiFi reads
- ONT (Ultra long)

Best case: both

Other data types can improve phasing and contiguity:

- Parental data
- Hi-C
- Strand-Seq
- Optical genome mapping

Technical background

Phased assembly of diploid genomes


Some commonly used assembly tools for human genomes:

Verkko [Rautainen et al. 2023], **Hifiasm** [Cheng et al. 2021], **Flye** [Kolmogorov et al. 2020]

Phased assembly of diploid genomes


Some commonly used assembly tools for human genomes:

Verkko [Rautainen et al. 2023], **Hifiasm** [Cheng et al. 2021], **Flye** [Kolmogorov et al. 2020]




Hifiasm: assembly with overlap graphs
[adapted from Li & Durbin 2024]

From unitigs to phased assembly




Hands on: let's assemble something

An example, HiFi-only assembly



An example, HiFi-only assembly



```
samtools view sample.bam "chr6:117000000-120000000" -b > reads_region.bam
```

```
samtools fastq sample_local.bam > reads_region.fastq
```

```
hifiasm reads_region.fastq -o my_assembly #<- here, we could include ONT, parental reads or Hi-C
```

(all code on github!) ²¹

gfa: graphical fragment assembly files

hifiasm outputs:

```
"main" .gfa output files
my_assembly.bp.hap1.p_ctg.gfa
my_assembly.bp.hap2.p_ctg.gfa
my_assembly.bp.p_ctg.gfa
my_assembly.bp.p_utg.gfa
my_assembly.bp.r_utg.gfa

other outputs:
my_assembly.bp.hap1.p_ctg.lowQ.bed
my_assembly.bp.hap1.p_ctg.noseq.gfa
my_assembly.bp.hap2.p_ctg.lowQ.bed
my_assembly.bp.hap2.p_ctg.noseq.gfa
my_assembly.bp.p_ctg.lowQ.bed
my_assembly.bp.p_ctg.noseq.gfa
my_assembly.bp.p_utg.lowQ.bed
my_assembly.bp.p_utg.noseq.gfa
my_assembly.bp.r_utg.lowQ.bed
my_assembly.bp.r_utg.noseq.gfa
my_assembly.ec.bin
my_assembly.ovlp.reverse.bin
my_assembly.ovlp.source.bin
```



```
S      utg000001l  TCCTGCACACCTATTTAAATGGCTAATACTAAAAAGATGACCTAGGAATTATTGGAGAGAATGTGAAACAATTGGAAACCTCATACATTGCTGGTG
S      utg000002l  AATACTATTACGCCATAAAAAGAAATGATATCATATTCTTTGCAGCAACAAGATGGAACTGGAGGCCATTATCCTAAGCTAACTAACACCGGGGCAG
S      utg000003l  CTGTTGTGCTAGCAAACTACTAGGCTTATTCACTTCTCTCAACTATTTTGTACCCATTAAACATCCATCACTTCTCTCCCACTAATCCCCACT
S      utg000004l  ACGATATTGATTCTTCTACCCATGAGCATGGAATGTTCTTCCATTGTGTTGATCTCTCTTTATTTCATTGAGCAGGGGTTGTAGTTCTCTTGA
S      utg000005l  AGAAAAATCAGCCCGTCCATCTTCTCTGATTATTTCCAGGGAAGACCTAGGCTCACATTAGCAGGGCTTAGCTGCTTCCGAGAGCAACAGGCGA
S      utg000006l  TGGTAGATTTTAGAATAAGTACCATTGTGGCACTCAGAAGAATGTATTTCTGTTGATTGGGCTAGAAAGTTCTGTAGACATCTACTAGGTCACCTT
S      utg000007l  GAGGAACTGGTACCATTCTTCTGAAACTATTTCCAATCAATAGAAAAAGAGGGAATCTCTCCTAACTCACTTTATGAGCCAGCATCACTTCTGTATC
S      utg000008l  GTTCCCATCATTTTGTAGGCCTCAGGATTAAATATGAAACCAGCCAAGAATACTAAACATTTGAGAAAGATTCTCACATTGAATCAGATTCCAAAA
S      utg000009l  GGGAACTCAGAATCTTCTGTAGGTGAATAACTGATATCTAAATTTAATGTTTGGGGAAAAGATATTTAAAAAAGATACTAGCCATATCATTACAT
S      utg000010l  TGTCAATTTGGGTGTTAAAAATAAGTCAAGCGGACTTAAAACTCTTACCCTACCAGAGAGAAAATTATTTAGAGCTACCTCATCAGATGTGCCTC
[...]
L      utg000001l  +      utg000003l  +      19822M  L1:i:85183      L2:i:0
L      utg000002l  +      utg000003l  +      15094M  L1:i:82821      L2:i:0
L      utg000003l  +      utg000004l  -      20551M  L1:i:153718     L2:i:0
L      utg000003l  +      utg000005l  -      9994M   L1:i:164275     L2:i:0
L      utg000003l  -      utg000001l  -      19822M  L1:i:154447     L2:i:0
L      utg000003l  -      utg000002l  -      15094M  L1:i:159175     L2:i:0
L      utg000004l  +      utg000003l  -      20551M  L1:i:9079       L2:i:0
L      utg000004l  -      utg000006l  -      16976M  L1:i:12654      L2:i:0
L      utg000004l  -      utg000007l  +      983M    L1:i:28647      L2:i:0
L      utg000005l  +      utg000003l  -      9994M   L1:i:26037      L2:i:0
L      utg000005l  -      utg000007l  +      17941M  L1:i:18090      L2:i:0
L      utg000006l  +      utg000004l  +      16976M  L1:i:83417      L2:i:0
L      utg000006l  -      utg000008l  +      17697M  L1:i:82696      L2:i:0
L      utg000007l  +      utg000008l  +      11045M  L1:i:66696      L2:i:0
L      utg000007l  -      utg000005l  +      17941M  L1:i:59800      L2:i:0
L      utg000007l  -      utg000004l  +      983M    L1:i:76758      L2:i:0
L      utg000008l  +      utg000009l  -      17683M  L1:i:1300       L2:i:0
[...]
```

S: Segment / Sequence

L: Link

Assembly graph visualization with *Bandage*

bandage -> load graph.... -> my_assembly.p_utig.gfa




nodes: Segments (“utigs”)
connections: Links



Excursion: interpreting assembly (unitig / contig) graphs

Excursion: Understanding assembly graphs (.gfa)

HiFi reads; 5Mbp repeat-rich region




Adverse signs:


- branching / tangles 
- 'tips': loose ends 

Excursion: Understanding assembly graphs (.gfa)



HiFi reads; 5Mbp repeat-rich region



HiFi reads; (peri) centromeric region on chr15




Adverse signs:


- branching / tangles 
- 'tips': loose ends 

Excursion: Understanding assembly graphs (.gfa)



HiFi reads; 5Mbp repeat-rich region




HiFi reads; (peri) centromeric region on chr15



Adverse signs:

- branching / tangles 
- 'tips': loose ends 



Excursion: Verkko assembly with $>100\times$ Hifi + $\sim 80\times$ ONT UL



Excursion over

Assembly graph visualization with *Bandage*



bandage -> load graph.... -> my_assembly.p_utig.gfa



nodes: Segments (“utigs”)
connections: Links

Assembly graph visualization with *Bandage*

bandage -> load graph.... -> my_assembly.p_utig.gfa




nodes: Segments (“utigs”)
connections: Links

no tips or tangles


manipulating gfa files - mapping segments back to the reference

```
gfatools gfa2fa unitigs.gfa > unitigs.fa
```

```
minimap2 -x asm5 -a hg38.fa unitigs.fa > aligned_utigs.sam
```




From unitigs to contigs: How to connect the sequences?




Additional data types (Parent info, Hi-C, ONT UL) will greatly improve this step

From unitigs to contigs: How to connect the sequences?




Additional data types (Parent info, Hi-C, ONT UL) will greatly improve this step

From unitigs to contigs: How to connect the sequences?



Additional data types (Parent info, Hi-C, ONT UL) will greatly improve this step


From unitigs to contigs: How to connect the sequences?



Additional data types (Parent info, Hi-C, ONT UL) will greatly improve this step


... and this is what `asm_hap1` and `asm_hap2` in our case are

```
gfatools gfa2fa unitigs.gfa > unitigs.fa  
minimap2 -x asm5 -a hg38.fa unitigs.fa > aligned_utigs.sam
```



- (1) Contiguity and Correctness with *Quast*
- (2) k-mer based evaluation with *Merqury*
- (3) Completeness with *BUSCO*

- (1) Contiguity and Correctness with *Quast*
- (2) k-mer based evaluation with *Mercury*
- (3) Completeness with *BUSCO*



Report

	bothhaps
# contigs (≥ 0 bp)	2
# contigs (≥ 1000 bp)	2
# contigs (≥ 5000 bp)	2
# contigs (≥ 10000 bp)	2
# contigs (≥ 25000 bp)	2
# contigs (≥ 50000 bp)	2
Total length (≥ 0 bp)	5939444
Total length (≥ 1000 bp)	5939444
Total length (≥ 5000 bp)	5939444
Total length (≥ 10000 bp)	5939444
Total length (≥ 25000 bp)	5939444
Total length (≥ 50000 bp)	5939444
# contigs	2
Largest contig	3044186
Total length	5939444
GC (%)	38.20
N50	3044186
N75	2895258
L50	1
L75	2
# N's per 100 kbp	0.00

```
quast.py <(cat asm_hap1.fa asm_hap2.fa)
```

- (1) Contiguity and Correctness with *Quast*
- (2) k-mer based evaluation with *Merqury*
- (3) Completeness with *BUSCO*

```
$ cat out.qv
label      erroneous-kmers  total-kmers  QV      Error-rate
asm_hap1   15                   3044166      66.296   2.34641e-07
asm_hap2   34                   2895238      62.5242  5.59214e-07
Both       49                   5939404      64.0576  3.92858e-07
```


```
meryl count k=21 ../../pulled_reads.fastq.gz output meryl_read_counts
$MERQURY/merqury.sh meryl_read_counts asm_hap1.fa asm_hap2.fa
```


Assembly assessment

- (1) Contiguity and Correctness with *Quast*
- (2) k-mer based evaluation with *Merqury*
- (3) Completeness with *BUSCO*

```
$ cat out.qv
label      erroneous-kmers  total-kmers  QV      Error-rate
asm_hap1   15                3044166     66.296  2.34641e-07
asm_hap2   34                2895238     62.5242 5.59214e-07
Both       49                5939404     64.0576 3.92858e-07
```

“k-mers which appear twice on our assembly:
usually appear ~30X in the raw reads”
-> homozygous peak



```
meryl count k=21 ../.../pulled_reads.fastq.gz output meryl_read_counts
$MERQURY/merqury.sh meryl_read_counts asm_hap1.fa asm_hap2.fa
```

- (1) Contiguity and Correctness with *Quast*
- (2) k-mer based evaluation with *Merqury*
- (3) Gene completeness with *BUSCO* (or its reimplementations 'compleasm')

BUSCO analysis (4.1.4)


Complete 99.1% (S+D)

● Single-copy 98.4%


● Duplicated 0.7%

● Fragmented 0.2%


○ Missing 0.7%




Complex regions require new visualizations




Complex regions require new visualizations




Miropeats plot [SVbyEye; [Porubsky et al. 2025](#)]




Complex regions require new visualizations




Miropeats plot [SVbyEye; Porubsky et al. 2025]




Dotplot [NAHRWhals; Höps et al. 2023]




Complex regions require new visualizations



Miropeats plot [SVbyEye; Porubsky et al. 2025]




Dotplot [NAHRWhals; Höps et al. 2023]




SV calling tools: e.g. *Dipcall*, *Smartie-sv*, *SVIM-asm*, *PAV*. [See e.g. review: Liu et al. 2024]

Example: Long duplications can obscure breakpoints




Deletion is visible, breakpoints are not.

Example: Long duplications can obscure breakpoints



assembly_h1 and _h2: mapped 'back' to hg38 to reveal more detail


The most complex genomic loci




[Logsdon et al. 2024]

Pan-genome graphs to understand complex variation


RDH;RHCE
GRCh38



[Liao et al. 2023]



Assemblies are only as good as the data - and always require expertise, time and careful QC.



Assemblies are only as good as the data - and always require expertise, time and careful QC.

It's worth it: De-novo assembly can answer questions that standard read-mapping can't.

Genome assembly algorithms

[1] Heng Li & Richard Durbin 2025: Genome assembly in the telomere-to-telomere era; Nat. reviews genetics

[2] Rautiainen et al. 2023: Telomere-to-telomere assembly of diploid chromosomes with Verkko

[3] Kolmogorov et al. 2019: Assembly of long, error-prone reads using repeat graphs

Assembly-based variant calling

[4] Olson et al. 2023: Variant calling and benchmarking in an era of complete human genome sequences

[5] Ebler et al. 2022: Pangenome-based genome inference allows efficient and accurate genotyping across a wide spectrum of variant classes

[6] Liu et al. 2024: Tradeoffs in alignment and assembly-based methods for structural variant detection with long-read sequencing data

Tutorials

Assembly step by step (highly recommend!):

https://training.galaxyproject.org/training-material/topics/assembly/tutorials/vgp_genome_assembly/tutorial.html

Assembly QC: (Quast, BUSCO, Mercury, Chromeister):

<https://training.galaxyproject.org/training-material/topics/assembly/tutorials/assembly-quality-control/tutorial.html>



https://github.com/WHops/ESHG2025_assembly_workshop